

ALTO 200 MOBILE FLOOR HOIST

The ALTO 200 is designed as a medium sized heavy duty mobile patient lifting hoist with a safe working load of 200kg with a compact ergonomic design. The ALTO provides superior patient lifting through the use of space link technology whilst remaining compact. The ALTO is the lifter of choice for institutional use or home use where a high capacity lift is required.

Features:

- Castors are swivel style, 100mm (4") diameter castors on the front and rear, with foot operated step-on lockable brakes on the rear
- Handheld control on an extendable cord with separate push buttons for raising and lowering the hoist at a single speed
- Battery mounted on the mast of the hoist, an alarm indicates low battery charge
- Battery charging is achieved by connecting the lifter to an off-board battery charger and then to a power source
- Emergency stop button
- Emergency lowering device
- **Space Link-** Increases maximum lifting height, enhances floor lift capability

Specifications:

- SWL 200kg
- Base Length: 1150mm
- External Base Width: 645mm
- Internal Base Width Closed: 578mm
- Internal Base Width Open: 1130mm
- Base Width Adjustment: Electric
- Under Bed Clearance: 125mm
- Minimum Height Spreader Bar: 362mm
- Maximum Height Lift Arm: 1785mm
- Machine Weight: 43kg
- Minimum Storage Height: 1350mm
- Average lift range: 1060mm
- Under base height: 50mm
- Maximum Hoist Reach: 875mm
- Minimum Hoist Reach: 445mm

Options:

- Tilt/Pivot frame
- Weighing unit
- 4 point yoke

