

CHAMELEON Collapsible Shower & Commode Chair

The CHAMELEON is the first Shower and Commode Chair developed specifically for travelling and comes in its own carry case. With its maximum width of 54cm the chair easily passes through standard-width toilet and shower doors, and adjusts to different height pedestals.

Features:

- Lightweight collapsible aluminium frame with interchangeable padded support seat and moulded toilet seat
- Collapses into carry bag, within international airline cabin baggage limits
- Four lockable castors (standard)
- Adjustable arm rests
- Dismantles and re-assembles **without tools**
- Flexible usage from private room to public areas

Optional Accessories:

- Two 18" pneumatic rear wheels with saddle-bag case
- Brakes for 18" wheels
- Foot rest for extra support
- Anti Pressure pillow
- Pan
- Backpack straps

Specifications:

- Manufactured in accordance with European standards
- **Packed in Suitcase:** Width:50cm, Depth:49cm, Height:17cm
- **Weight packed in suitcase:** 9.9kg
- Seat Width: 46cm
- Seat Depth: 42cm
- Height Range: 45cm – 55 cm
- User Weight: 120kg

