

FLEX-LIFT Collapsible Mobile Floor Hoist

The **FLEX-LIFT™** is a mobile floor hoist designed for travelling. It can be folded away in seconds, fits into a regular car boot.

Features:

- Lightweight metal frame
- Picks up from the floor with a lifting range of 116cm
- Frame folds for easy storing when not in use
- Assembled and disassembled without the use of tools
- Can be easily carried by one person
- Lockable castors
- Easy carried in 2 bags – 13kg & 14kg

Specifications:

Max Lifting weight: 130kg
Total weight: 23.9kg

Flex Lift assembled:

Height: 117-174cm
Hoisting range: 116cm
Yoke, lowest point: 46cm
Yoke, highest point: 162cm
Width, folded: 28cm
Width, outside: 77.5cm
Width, inside: 66cm
Clearance, undercarriage: 5.5cm

Undercarriage:

Length: 113cm
Width folded: 22.5cm
Weight in bag: 14.3kg

Hoisting Mast:

Hoisting Mast length: 110cm
Weight in bag: 13.9kg

Battery:

24 volt (25-30 lifting cycles)
Charges through hand control at external charging station
Charging time 4-5 hours

