

## GH2F Portable Ceiling Hoist System

**GULDMANN Ceiling Hoists** offer full flexibility for the lifting, transferring and transport of elderly and disabled persons in private homes, hospitals, institutions and swimming pools.

Additionally the GULDMANN lifting module comes with the option of being fully portable. This portable option makes GULDMANN easy to move from one room to another as it installs and removes itself from the rail system automatically. With the help of the specially designed transport and storage trolley, it can be moved without putting a strain on staff.

The **GULDMANN** hoist is suitable for mounting in single-track, multi-track room covering and combined rails systems. A full range of lifting slings are available and the hoist system can handle all lifting situations.

By combining this hoisting system with PARA-MOBILITY's expertise in designing and manufacturing engineering structures, we can accommodate the specific needs required by clients.

### Features include:

- High lifting capacity - 200kg to 250kg
- Quiet operation - will not disturb others or attract attention
- Battery operated with dual, powerful, long lasting batteries - has up to 140 lifts before automatic recharging occurs
- Emergency lowering easily activated by the user or carer
- Designed for use in damp environments
- A wide range of track combinations - straight, curved and switching, which allows for extensive design flexibility.

### Specifications:

- Lift range 2 metres
- Hand control operation
- Manufactured in accordance with ISO 9999 : 12 36 12.
- Tracking materials: lightweight aluminium; max span 8 metres
- GH2F version is fully portable and removes itself from the rail system automatically.



View a video of GH2F hoist in action: [www.paramobility.com.au/products/track\\_ceiling.asp](http://www.paramobility.com.au/products/track_ceiling.asp)

### Manufacturer and Importer of Specialised Disability Equipment

**Head Office:** Ph: (02) 9653 2719 Fax: 02 9653-3719 Email: [sales@paramobility.com.au](mailto:sales@paramobility.com.au)

**QLD Office:** Ph: (07) 5549 3106 Fax: (07) 5549 1313 Email: [qld@paramobility.com.au](mailto:qld@paramobility.com.au)

[www.paramobility.com.au](http://www.paramobility.com.au)

A division of Jaswil Pty Ltd – ABN 29 081 258 291