

The IBIS Vehicle Access Lifter system has been designed and manufactured by Para Mobility to allow a disabled person to be seated in the front passenger seat of most cars, 4WD vehicles, the side entry of vans and caravans. It is operated by the person or carer, with a press button hand control allowing the person to be transferred from a wheel chair into the vehicle **WITHOUT ANY LIFTING**, whilst being completely portable.

Features:


- Operation can be by either the user or the carer
- Mounting bracket is neatly located beside the front door column
- Press button hand control which is powered from the vehicle battery
- Can be easily removed and stored in the boot when not in use and whilst travelling
- Lightweight for easy handling - 11kgs
- Includes Ibis sling, tailor made slings can be arranged

Specifications:

- Manufactured by Para Mobility in accordance with ASA 3581
- Lifting capacity 100kgs
- Mechanical Actuator Lift mechanism

Optional Accessory:

- Carry/storage bag



View video of the Ibis in action: www.paramobility.com.au/products/ibis_vehicle.asp

Manufacturer and Importer of Specialised Disability Equipment

Head Office: Ph: (02) 9653 2719 Fax: 02 9653-3719 Email: sales@paramobility.com.au

QLD Office: Ph: (07) 5549 3106 Fax: (07) 5549 1313 Email: qld@paramobility.com.au

www.paramobility.com.au

A division of Jaswil Pty Ltd – ABN 29 081 258 291