

Mobile Height Adjustable Change Table

Designed and manufactured by Para Mobility (Australian owned and operated), the mobile height adjustable change table allows easy transfer from a wheel chair when the table top is at its lowest point. The table can then be raised to a working height for the carer. The table top is removable and can be tailor made to suit the dimensions of the room and the user.

NB: Care should be taken in a wet environment.

Features:

- Lowers to allow a person to slide transfer or use a floor lifter from a wheelchair to the table top
- Hand held control for easy operation
- Battery / electric operation. Multiple lifts between charging
- Restraint strap included
- Table top can be easily removed for easy storage, transport or fitting into tight areas
- Padded and covered with high grade PVC material for urine and stain control
- Unit is on castors with 2 locking wheels
- Framework is electroplated (zinc plating) and powder coated for corrosion resistance
- Used for washing and changing in either a wet or dry area
- Designed and built in accordance with ASA 3581

Specifications:

- Standard table top size: 1850mm long x 700mm wide. Sizes can be made to suit the user
- Rises to approximately 930mm
- Lowest point – approximately 440mm
- Drop down sides: 150mm standard high x 1200 long. Sides can be made to order
- Has clearance of 135mm at bottom allowing floor hoist to fit underneath
- SWL (safe working load) 150kg

Optional Extras:

- Drop down sides
- Additional straps
- Foot Pedal operation



Manufacturer and Importer of Specialised Disability Equipment

Head Office: Ph: (02) 9653 2719 Fax: 02 9653-3719 Email: sales@paramobility.com.au

QLD Office: Ph: (07) 5549 3106 Fax: (07) 5549 1313 Email: qld@paramobility.com.au

www.paramobility.com.au

A division of Jaswil Pty Ltd – ABN 29 081 258 291