

PLATYPUS Beach/All Terrain Wheelchair

The Platypus Beach / All Terrain Wheel Chair has been designed to allow access to areas where conventional wheelchairs are not able to go, such as beaches, walking trails and many backyards. The Platypus Beach / All Terrain Wheel Chair clears standard doors and gates increasing its accessibility. The Platypus Beach and All Terrain Wheel Chair is suitable for people with disabilities of all sizes from children to adults. Naturally the heavier the person the more energy will be required for pushing. The Wheeleez wheels are attached with a quick release mechanism eliminating lost axle pins. When disassembled no component weighs more than 13kg, making it easier to handle and pack away for transport.

Features:

- Lightweight
- Easy to manoeuvre
- Swing away arms for easy transfer
- Can negotiate a standard door or gate
- Corrosion resistant
- Open weave fabric to allow cool air flow
- SWL 120kg
- Specially designed Wheeleez Wheels

Specifications:

- Seat Height 57cm
- Seat Width 44cm
- Chair Height 100cm
- Chair Width 68cm
- Chair Length 97cm
- Weight when assembled 35kg

Optional Extras:

- Child adaptor
- Choice of wheels
- Lap Belt

