

## PLATYPUS Mobile Change Table

Platypus Products Australia has designed this change table as a lightweight alternative to standard change tables. The unique design allows a person to be changed in any area (wet or dry) whilst being completely mobile.

### Features:

- Sturdy frame
- 700mm (wide) x 1700mm (long) x 800mm (high)
- Stainless steel castors (2 lockable)
- SWL 80kg
- Accidental burns are prevented as plastic does not absorb the heat
- Open weave fabric which allows the easy passage of air and water
- Easy to dry and clean
- Rounded surfaces which improve comfort and reduces the risk of skin tears

### Optional Accessories:

- Safety Belt/s
- Raised Head
- Storage area

