

We manufacture a wide range of light weight commode chairs, shower chairs and bedside commodes. The Platypus range of chairs are made from furniture grade plastic providing a very strong, durable and lightweight product comfortable for the client and easy to use for the carer.

Features:

- Lightweight
- Cush n Soft Padded seat (open, closed or junior size)
- Rounded surfaces which are comfortable to hold and helps to reduce skin tears
- Accidental burns are prevented as plastic does not absorb heat
- Open weave fabric allows easy passage of air and water
- Easy to dry and clean
- Chairs are easily manoeuvred with the combination of quality castors and a lightweight frame
- Swing away arms for side transfer
- Folding height adjustable footrests
- SWL 120kg

Optional Accessories:

- Leg rest (single or double)
- Bariatric version
- Child adaptor
- Safety arms
- Pan holder (holds standard bed pan and easily clicks on)
- Safety belts (harness or lap)
- Self propelling chair with 24" wheels (includes user operated brakes)



Standard Shower/Commode Chair with open Cush n Soft seat, harness belt and pan holder



Standard Shower/Commode chair with single leg rest and child adaptor



Standard Shower/Commode chair with double leg rest



Self propelling Shower/Commode Chair with safety arms