

SALSA 200 MOBILE STAND UP HOIST

The SALSA 200 is designed as an active standing aid allowing a single carer to raise a patient into a comfortable secure standing position. The SALSA is ideal for transferring from bed to chair, or chair to commode. Allied to suit the SALSA, is a padded sling which allows superior access to remove clothing for toileting purposes. The SALSA is equally at home in an institutional or home environment.

Features:

- Castors are swivel style, 100mm (4") diameter castors on the front and rear, with foot operated step-on lockable brakes on the rear
- Handheld control on an extendable cord with separate push buttons for raising and lowering the hoist at a single speed, and opening and closing the base legs. There is a hook on the handpiece for mounting it to the hoist when not in use
- Battery mounted on the mast of the hoist, an alarm indicates low battery charge
- Battery charging is achieved by connecting the lifter to an off-board battery charger and then to a power source
- Removable footrest with a slip resistant surface and a raised edge
- Knee pad is padded, upholstered, contoured, height adjustable and comes with a neoprene calf support strap with four points of key-hole attachments
- Emergency stop button
- Emergency lowering device

Specifications:

- SWL 200kg
- Base Length: 1050mm
- External Base Width: 655mm
- Internal Base Width Closed: 585mm
- Internal Base Width Open: 1035mm
- Base Width Adjustment: Electric
- Under Bed Clearance: 125mm
- Minimum Height Lift Arm: 825mm
- Maximum Height Lift Arm: 1670mm
- Machine Weight: 48.2kg
- Minimum Storage Height: 1200mm

Options:

- Velcro or buckle chest/stomach strap on sling (come in 2 sizes)



Manufacturer and Importer of Specialised Disability Equipment

Head Office: Ph: (02) 9653 2719 Fax: 02 9653-3719 Email: sales@paramobility.com.au

QLD Office: Ph: (07) 5549 3106 Fax: (07) 5549 1313 Email: qld@paramobility.com.au

www.paramobility.com.au

A division of Jaswil Pty Ltd – ABN 29 081 258 291