

## SONATA 150 MOBILE FLOOR HOIST

A compact mobile battery operated hoist which can be used to lift and transfer a person to and from a bed, chair and the floor. The hoist is compatible with a standard spreader bar or tilting /pivot frame. Leg opening can be operated both manually and electronically. Suitable for both domestic use and healthcare facilities.

### Features:

- Epoxy coated steel frame
- Upright mast with an angled boom and a removable padded spreader bar
- The standard spreader bar has two points of attachment for a sling
- The hoist features a "Space Link" device, which increases maximum lifting height
- The U-shaped base has manual leg adjustment via a large lever with a contoured handle
- 3 positions for legs
- Dismantles into two pieces
- Hand held on an extendable cord with separate push buttons for raising and lowering the hoist at a single speed. There is a hook on the hand piece to enable it to be located on the hoist when not in use
- 12 month warranty

### Specifications:

- SWL 150kg
- The hoist has an on-board battery charger. This requires the hoist to be stored with access to a power point
- Emergency stop button

### Measurements:

Minimum height of spreader bar	340mm
Maximum height of spreader bar	1685mm
External base width (closed)	605mm
Internal width of base (closed)	530mm
Internal width of base (open)	1010mm
Base height (clearance required)	125mm
Length of base	1005mm
Overall weight	39kg

### Options:

- Tilting/Pivot frame
- A range of slings can be selected to use with the hoist in sizes ranging from extra small to XXL.
- Slings available include: general purpose, general purpose with head support, tilt/pivot, amputee and paediatric slings



### Manufacturer and Importer of Specialised Disability Equipment

**Head Office:** Ph: (02) 9653 2719 Fax: 02 9653-3719 Email: sales@paramobility.com.au

**QLD Office:** Ph: (07) 5549 3106 Fax: (07) 5549 1313 Email: qld@paramobility.com.au

www.paramobility.com.au

A division of Jaswil Pty Ltd – ABN 29 081 258 291