

Designed and manufactured by Para Mobility the therapy table is electronically operated, has press button height adjustable controls and can be made to any size required.

Features:

- Lowers to allow a person to slide transfer or use a lifter from a wheelchair to the table top
- Hand held control for easy operation
- Battery/electric operation with easy overnight recharging
- Table top is removable for easy storage
- Padded and covered with high grade PVC for urine and stain control
- Mobile with 4 locking castors
- Framework is electroplated and powder coated to a colour that suits
- Used for washing, therapy and changing in either wet or dry areas
- Manufactured to any size as required
- Designed and built in accordance with ASA 3581

Optional Extras:

- Drop Down Sides
- Restraint straps



Manufacturer and Importer of Specialised Disability Equipment

Head Office: Ph: (02) 9653 2719 Fax: 02 9653-3719 Email: sales@paramobility.com.au

QLD Office: Ph: (07) 5549 3106 Fax: (07) 5549 1313 Email: qld@paramobility.com.au

www.paramobility.com.au

A division of Jaswil Pty Ltd – ABN 29 081 258 291