

Tilt Tables/ Standing Frames

Designed and manufactured by Para Mobility with input from Physiotherapists and Occupational Therapists for use in special schools, hospitals and home use. The tilt tables may be powered direct to a 240v GPO, or by rechargeable battery.

Features:

- Can be table tailor-made to suit size or growth of the user
- Hand held pendant control
- Four locking castors
- Electroplated and powder coated colour to suit
- Child size with an adult extension available.
- Easy transfer direct from wheelchair or lifter to the horizontal table
- Safe Working Load of 150kg
- Size: 1560mm long x 650wide (suit Child) Head rest can be added 400mm long to suit adult making length 1900mm x 650wide
- Adjustable, wide and comfortable, chest / pelvic and split leg support straps
- Totally adjustable padded side supports for knees, legs, hips and upper chest to suit the users requirements.(for Work Station Model)

Optional Extras:

- Adjustable ankle, side, chest, and knee supports
- Additional straps can be added



Basic Model –Adult size



Work Station Model