

Wall Mounted Height Adjustable Change Table

Designed and manufactured by Para Mobility the Wall Mounted Height Adjustable Change Table can be used in schools, private homes, community centres, swimming pool change rooms and leisure centres.


The Wall Mounted Height Adjustable Change Table has been designed to be stored flat against the wall within its own frame when not in use, creating extra space. When ready for use, the table is easily lowered to the desired height with no power required for operation.

Features:

- Lightweight aluminium construction
- Velcro secured cushions ensure the user's safety
- Table top and surrounding cushions are covered in Medi-Flex for stain minimisation
- Three (3) height positions with easy adjustment. 75mm apart
- Fitted with one (1) restraint strap for safety. Additional strapping can be added if necessary

Specifications:

- Manufactured by Para Mobility in accordance with ASA 3581 standards, with lightweight aluminium frame
- Three (3) height adjustments – 75mm apart
- Can be mounted on brick or stud walls, full installation instructions are supplied on delivery
- **Recommended** installation height: 600mm from floor to frame
- **Standard bed size:** 1850mm long x 750mm wide; sizes can be tailor made to user requirements. **SWL - 120kg**
- **Small bed size:** 1200m long x 750mm wide. **SWL - 80kg**
- NB - an additional 120mm needs to be added to accommodate the frame when checking available space



Manufacturer and Importer of Specialised Disability Equipment

Head Office: Ph: (02) 9653 2719 Fax: 02 9653-3719 Email: sales@paramobility.com.au

QLD Office: Ph: (07) 5549 3106 Fax: (07) 5549 1313 Email: qld@paramobility.com.au

www.paramobility.com.au

A division of Jaswil Pty Ltd – ABN 29 081 258 291