

Kingfisher Pool Hoist by Para Mobility

Introducing the newest member of our collection.... The Kingfisher Pool Hoist - Manufactured in Australia by Para Mobility and released in Feb 2014 the Kingfisher Pool Hoist will allow access to a pool where traditional methods are not able to be used. In close consultation with Occupational Therapists and other industry professionals the Kingfisher Pool Hoist is a simple and safe operation that is both comfortable and dignified. The Kingfisher is ideal for home pools, retirement villages and hotels etc. It is user-operated and secured to the pool deck with minimal fuss. A trolley is also available which allows it to be easily moved between locations and/or storage. Manufactured to Australian Standard — AS3581

Key Features

The Kingfisher Pool Hoist has a Safe Working Load (SWL) of 150kg, rechargeable Linak 24v battery, Floating hand control, Load bearing Lift up arm rests for ease of transfer, removable and adjustable foot rest, flip up seat for storage, wheel away transport trolley option for storage, lap belt for safety, left or right side mounting, extra heavy duty Linak Actuator with emergency lifting.

Specifications

• Lifting Capacity:	SWL 150kg
• Weight of Hoist:	48kg (without the battery)
• Max—Deck to water level:	up to 250mm
• Max—Set Back:	500mm
• Stow away size:	H 1230 x W 550 x D 400
• Seat size:	400mm x 400mm
• Back rest size:	400mm 300mm
• Foot rest size:	300mm x 150mm (can be modified)
• Powder coated:	White

Materials

• Base Plate:	316g Stainless Steel
• Lifting Frame:	316g & 304g Stainless Steel
• Actuator:	Linak
• Battery:	Linak
• Control Box:	Linak
• Floating Hand Control:	Linak
• Seat:	High Density Polyethylene
• Seat Belt:	Interwoven nylon/polyester
• Foot rest:	High Density Polyethylene
• Hand Control:	Floating –Linak

Standard Hoist Kit inclusions

- Hoist
- Lift up arm rests
- 24v Battery / Battery charger / cover
- Hand control
- Seat belt
- All weather cover
- Mounting Kit

Optional

- Foot rest, Chest strap, head rest, Transport trolley

