

Waverley Glen C450/C625 Ceiling Hoist System

The Waverley Glen C Series Ceiling Hoist system delivers innovation, affordability and quality. The C series has been designed with a (delete) many standard features which allow it to be operated by one carer in many environments including: private homes, hospitals, institutions and swimming pools.

The units are very compact yet powerful enough to lift up to 280kg, the C series will be an asset to any home or care environment.

Features include:

- Lifting capacity 200kg (C450) or 280kg (C625)
- Quiet operation
- Quick lifting speed to reduce patient transfer times
- Emergency lowering easily activated by the user
- Emergency stop pull cord
- Audible low battery indicator and visual battery/charge level display
- Digital display that indicates number of lifts, battery levels and lift status
- Auto shut off when not in use to conserve battery life
- Manual emergency lowering/lifting device
- Pneumatic hand control which is completely waterproof
- Single Rail, Curves, room covering systems & turntables

Specifications:

- Lift range 213cm, allows hoisting from the floor
- Hand control operation
- Lift Motor- 24v DC
- Traverse Motor- 24v DC – applicable in power traverse models only
- Batteries – 2x12v DC
- Unit weight under 10kg

Optional Accessories:

- Power traverse with 'return to charge'
- Automatic return to charge
- Large range of slings



Manufacturer and Importer of Specialised Disability Equipment

Head Office: Ph: (02) 9653 2719 Fax: 02 9653-3719 Email: sales@paramobility.com.au

QLD Office: Ph: (07) 5549 3106 Fax: (07) 5549 1313 Email: qld@paramobility.com.au

www.paramobility.com.au

A division of Jaswil Pty Ltd – ABN 29 081 258 291