



# IBIS™ PORTABLE CAR ACCESS LIFTER



## Take control and enjoy greater independence

The IBIS™ Portable Car Access Lifter system has been designed and manufactured in Australia by Para Mobility to allow people with limited mobility the option of getting in a car, SUV, 4WD or van. Passengers can be seated in the front passenger seat of most cars, 4WD vehicles and the side entry of vans and caravans. It is operated by the person or carer, with a press button action allowing the person to be transferred from a wheelchair into the vehicle without any lifting, whilst being completely portable.

To further make the IBIS™ Access Lifter even more portable an IBIS™ Hinge Adapter Bracket (IHAB) is available, custom made to attach to the hinges of your car door, which eliminates the need for any permanent car modifications. Clip the IHAB onto the door hinge of the car, attach the IBIS™ to the bracket as normal, plug in the power source using the portable battery pack and you are ready to go.



# IBIS™ PORTABLE CAR ACCESS LIFTER **PRODUCT SNAPSHOT**

## FEATURES

- Operation can be by either the user or the carer
- Mounting bracket is neatly located beside the front door column
- Press button hand control which is powered from the portable battery pack or vehicle battery
- Can be easily removed and stored in the boot when not in use and whilst travelling
- Lightweight for easy handling - 11kgs
- Includes IBIS™ sling (tailor made slings can also be arranged)

## SPECIFICATIONS

- Manufactured by Para Mobility in accordance with ASA 3581
- Lifting capacity 100kgs
- Mechanical Actuator Lift mechanism

## INCLUSIONS

- Carry/storage bag
- Sling
- Battery pack
- IBIS™ hinge adapter bracket (IHAB)

No car modification required



View in action:

[www.paramobility.com.au/mobility-products/mobility-range/ibis-car-access-hoist](http://www.paramobility.com.au/mobility-products/mobility-range/ibis-car-access-hoist)

**Para Mobility**  
*Making everyday living easier since 1983*

Phone (02) 9651 4446

Delivery, installation + servicing Australia wide  
[sales@paramobility.com.au](mailto:sales@paramobility.com.au)

**[www.paramobility.com.au](http://www.paramobility.com.au)**

