



POWERED WALL MOUNTED HEIGHT ADJUSTABLE CHANGE TABLE



The power to get out and about!

Designed and manufactured in Australia by Para Mobility, the Power Wall Mounted Height Adjustable Change Table can be used in schools, private homes, community centres, swimming pool change rooms and leisure centres.

The Power Wall Mounted Height Adjustable Change Table has been designed to be stored flat against the wall within its own frame when not in use, creating extra space. When ready for use, the table is easily lowered to the desired height with the hand control – even low enough for a person to slide transfer from a wheelchair!



POWERED WALL MOUNTED HEIGHT ADJUSTABLE CHANGE TABLE PRODUCT SNAPSHOT

MEETS 'CHANGING PLACES' GUIDELINES

FEATURES

- Lightweight aluminium construction – powder coated white.
- Table top and surrounding cushions are covered in Medi-Flex for stain minimisation. Medi-Flex is a flame resistant fabric which is slip resistant, waterproof, dust mite, urine and stain resistant with an anti-bacterial additive.
- 700mm lift range.
- Height adjustable with or without a person on it.
- Twin Linak actuators and hand control.
- Fitted with one (1) restraint strap for safety. Additional strapping can be added if necessary.
- Manufactured by Para Mobility in accordance with ASA 3581 standards, with lightweight aluminium frame.
- Easily retro-fitted to any masonry or timber stud wall. Full installation instructions are supplied on delivery.

MEASUREMENTS


- **Standard bed size:** 1850mm long x 750mm wide; sizes can be tailor made to user requirements. SWL - 200kg.
- NB: Allow 125mm on either side of bed for frame when measuring available space.


Para Mobility
Making everyday living easier since 1983

Phone (02) 9651 4446

Delivery, installation + servicing Australia wide
sales@paramobility.com.au

www.paramobility.com.au



Registered NDIS Provider

