

PLATYPUS BEACH WHEELCHAIR



Everyone should enjoy the beach and its benefits

The Platypus Beach and All Terrain Wheelchair has been designed to allow access to areas where conventional wheelchairs are not able to go, such as beaches, walking trails and many backyards.

The Platypus Beach and All Terrain Wheelchair clears standard doors and gates increasing its accessibility. The Platypus Beach and All Terrain Wheelchair is suitable for people with disabilities of all sizes from children to adults.

The Wheellez wheels are attached with a quick release mechanism eliminating lost axle pins. When disassembled no component weighs more than 13kg, making it easy to handle and pack away for transport.



PLATYPUS BEACH WHEELCHAIR PRODUCT SNAPSHOT

FEATURES

- Lightweight
- Easy to manoeuvre
- Lap Belt
- Swing away arms for easy transfer
- Can negotiate a standard door or gate
- Corrosion resistant
- Open weave fabric to allow air and water flow
- SWL 120kg
- Specially designed Wheeliez Wheels

SPECIFICATIONS

- Seat Height 57cm
- Seat Width 44cm
- Chair Height 100cm
- Chair Width 68cm
- Chair Length 97cm
- Weight when assembled 35kg

OPTIONAL EXTRAS

- Child adaptor
- Choice of wheels
- Harness Belt

**Para Mobility**
Making everyday living easier since 1983

Phone (02) 9651 4446

Delivery, installation + servicing Australia wide
sales@paramobility.com.au

www.paramobility.com.au



Disability Discrimination Act 1992

Proud supporters of



Registered NDIS Provider

