



# LUNA CEILING HOIST



## Modern and dynamic, for greater flexibility

The Luna Ceiling Hoist with its unique and modern spherical shape and excellent functional features, is a true example of innovative Danish design. Luna is a fixed and portable ceiling hoist in one unit. Changing Luna from fixed to portable mode, and vice versa, can be done within seconds by simply detaching and reattaching the unit to the track. No tools needed. This unique two-in-one feature makes Luna the preferred ceiling hoist for care facilities as it provides a flexibility of choice.

Combine Luna with Ergolet's extensive program of ceiling and wall mounted tracks, and you have a ceiling hoist system providing complete solutions to meet the latest requirements in the care sector. Other features include Ergolet's patented spreader bar and sling hook system ensuring safe and simple sling attachment. The comfortable ergonomic hand control, with its clear function buttons, is designed to help make patient transfer a safe and simple task for carers.



# LUNA CEILING HOIST PRODUCT SNAPSHOT

## FUNCTIONS

- Class leading state-of-the-art design and functionality
- Weight capacity of 200kg and 275kg (patented Ergolet sling hooks - simple and safe to use)
- Fixed or portable in a modular design
- Charging through the hand control and wall/table charger
- Emergency lowering function and brake system
- Soft start and stop motor control
- Electrical emergency lowering and stop
- Low battery indicator (visual LED)
- Charging and power on indicators
- Overload protection



## SPECIFICATIONS

- Max. capacity: 200kg / 275kg
- Length of lifting strap: 2200mm
- Weight: 6.7kg w/o spreader bar
- Measurements: 320 x 130mm
- Battery operation: 36V lead acid batteries
- Charging: 100-230V. 0.5A
- Chassis: Aluminium
- Housing: High impact, moulded ABS
- Noise: Lwa during lifting/lowering without load is 49dB (A)  
Lwa during lifting/lowering at maximum load is 53dB (A)
- Lifting speed down: 44mm/ sec with 89kg load
- Lifting speed up: 30mm/ sec with 89kg load
- No. of lifts per charge: 50 lifts of 0.5m with 100kg
- ISO: ISO 10535
- Versions: Luna 200kg (white/grey) Luna 275kg (white/grey)  
Luna room to room 200kg (white) or (grey)

## ACCESSORIES

- Complete selection of 2, 3 and 4 point spreader bars to suit standard and bariatric needs
- Complete range of slings
- Weigh scales

**Para Mobility**  
*Making everyday living easier since 1983*

Phone (02) 9651 4446

Delivery, installation + servicing Australia wide  
[sales@paramobility.com.au](mailto:sales@paramobility.com.au)

[www.paramobility.com.au](http://www.paramobility.com.au)

